

One inflation device for all your dilation procedures.



Sphere[®] INFLATION DEVICE

Ease of use


The inflation device has a transparent barrel so that you can clearly see the contents and volume. The device has an ergonomic design. You can unlock the device with one hand, even when the device is at full pressure.

Controlled inflation and deflation

You can rapidly inflate or deflate your balloon, or you can flip a switch to engage the locking mechanism and make fine adjustments to the pressure.

Easy-to-read pressure gauge

The pressure gauge is visible in low light and displays pressure in atmospheric pressure (atm) and pound-force per square inch (psi). The gauge also indicates vacuum.



Sphere®

INFLATION DEVICE

During minimally invasive balloon dilation procedures, the Sphere Inflation Device can be used to inflate the balloon with fluid, to monitor the pressure in the balloon during the procedure, and to deflate the balloon after dilation.

Order Number	Reference Part Number	Fluid Capacity mL	Maximum Pressure atm
G34903	CID-60-15	60	15

Contact your local Everis representative or Customer Service for details.

Everis Medical

555 East Eliza St, Ste. A
Schoolcraft, MI 49087 USA
t: 888 902 2239
www.everismedical.com

Printed in U.S.A.
DOC-6552 Rev 1
Copyright © 2025 Everis

Everis™