


# Treat obstructive salivary gland disorders.



## SialoCath™

**SALIVARY DUCT CATHETER**

Illustration by Lisa Clark

### Scope-free irrigation


The catheter allows for irrigation without the use of a scope as a first-line, in-office treatment approach.

### Atraumatic design

The soft tip, flexible catheter provides a gentle treatment option for both in-office and surgical procedures.

### Salivary use

The catheter is specifically indicated for salivary duct irrigation.



# SialoCath™


## SALIVARY DUCT CATHETER

Used to aid in the removal and treatment of salivary duct stones and other duct conditions through flushing and/or fluid irrigation of the salivary ducts.

This flexible, single-lumen catheter can be inserted into the salivary duct either over a wire guide or through an appropriately sized introducer sheath.

Order Number	Reference Part Number	Catheter OD mm/inch	Catheter ID mm/inch	Compatible Wire Guide Diameter inch	Compatible Catheter Length cm	Tip Configuration
G35166	SDC-2.2-20	.74/.029	.53/.021	.018	15	Straight


## Related products



### Kolenda Salivary Access Introducer Set

Used for the introduction of diagnostic or interventional devices in the submandibular salivary duct.

Order Number	Reference Part Number	Introducer ID mm/inch	Wire Guide Diameter inch	Sheath Length cm
G24129	SDAS-6.0-18	2.2/.087	.018	5.5
G24132	SDAS-5.0-18	1.9/.074	.018	5.5



### Wire guides

Used to gain and maintain access to the salivary duct.

Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm	Tip Configuration
G02274	C-SF-15-50	.015	50	Straight

Some products or part numbers may not be available in all markets. Contact your local Everis sales representative or Customer Service for details.

#### Everis Medical

555 East Eliza St, Ste. A  
Schoolcraft, MI 49087 USA  
t: 888 902 2239  
www.everismedical.com

Printed in U.S.A.  
DOC 6574 Rev 1  
Copyright © 2025 Everis

**Everis™**