

Everis Sialendoscopy Compatibility Chart	Outer Diameter (mm)	5 Fr Kolenda Introducer* G24132	6 Fr Kolenda Introducer* G24129	0.015" Wire Guide G02274	0.018" Wire Guide G02255	NCircle® Basket G23797	NGage® Basket G23798	Advance® Salivary Duct Balloon Catheter Set G34991, G34992, G34993, G34994
STORZ ERLANGEN Miniature Endoscopes								
Miniature Straight Forward Telescope 0°	0.8	●	●					◐
	1.1	●	●	●				◐
	1.6	●	●	●	●	●	●	
STORZ MARCHAL Miniature Endoscopes								
Miniature Straight Forward Telescope 0°	0.75	●	●					◐
	0.89	●	●					◐
	1.1	●	●	●				◐
	1.3	●	●	●	●	●		◐
	1.6	●	●	●	●	●	●	
Examination Sheath	1.1	●	●					◐

Ensure fit/setup prior to use and following hospital procedures.

- = Compatible with scope
- ◐ = Compatible alongside a scope through the 6 Fr Kolenda Introducer Sheath
- * Compatible for scope through the Kolenda Introducer Sheath

© Everis 2025 DOC 6575 Rev 1