


Biodesign®

OTOLOGIC BUTTERFLY GRAFT



Harvest *results*,
not patient tissue.¹

A transformative approach to traditional
butterfly graft tympanic membrane repair^{1,2}

Reliable Closure

Biodesign® material remodels into natural host tissue with an overall success rate of 91% across published literature³⁻¹¹ with no statistically significant difference in audiometric results when compared to temporalis fascia.^{3, 12}

Excellent Handling

Biodesign® material is easy to manipulate, allowing for improved surgical precision during graft placement.³

Time Saving

The Biodesign® Otolgic Butterfly Graft is a new, off-the-shelf option for a minimally invasive procedural approach. It reduces the need to harvest autologous tissue, significantly decreasing intraoperative time.³



PROCEDURE
VIDEO

Everis™



Available product sizes
Shown at actual size.



Biodesign® Otologic Butterfly Graft

Order Number	Reference Part Number	Size cm	Nominal Thickness mm
G60285	ENT OTO BFLY 0.4 0.6	0.4, 0.6	0.25

Tips to help get the best possible results:

- Ensure adequate blood supply.
- Size the graft to allow some tissue overlap.
- Place the graft dry or hydrate it for less than one minute before placement.

References

1. Internal white paper; "Evaluation of industry opinion on endoscopic transcanal TM repair as a conservative and minimally invasive procedure". Available on request.

2. https://www.accessdata.fda.gov/cdrh_docs/pdf23/K232646.pdf; Published K232646 510k summary describing device usability testing data in benchtop and animal models following an endoscopic transcanal approach.

3. D'Eredità R. Porcine small intestinal submucosa (SIS) myringoplasty in children: a randomized controlled study. *Int J Pediatr Otorhinolaryngol.* 2015;79(7):1085-1089.

4. Cass ND, Hebbe AL, Meier MR, et al. Pediatric primary tympanoplasty outcomes with autologous and non-autologous grafts. *Otol Neurotol.* 2022;43(1):94-100.

5. Chen CK, Hsieh LC. Clinical outcome of exclusive endoscopic tympanoplasty with porcine small intestine submucosa in 72 patients. *Clin Otolaryngol.* 2020;45(6):938-943.

6. Barron C, Lukens J, Niermeyer W, et al. Investigation of novel grafts in use for pediatric tympanoplasty. *Ann Otol Rhinol Laryngol.* 2019;128(12):1111-1115.

7. Redaelli De Zinis LO, Berlucchi M, Nassif N. Double-handed endoscopic myringoplasty with a holding system in children: preliminary observations. *Int J Pediatr Otorhinolaryngol.* 2017;96:127-130.

8. James AL. Endoscope or microscope-guided pediatric tympanoplasty? Comparison of grafting technique and outcome. *Laryngoscope.* 2017;127(11):2659-2664.

9. Ranguis SC, Leonard CG, James AL. Prospective comparison of pediatric endoscopic lateral graft and interlay tympanoplasty. *Otol Neurotol.* 2021;42(6):867-875.

10. Wang N, Isaacson G. Collagen matrix as a replacement for Gelfilm for post-tympanostomy tube myringoplasty. *Int J Pediatr Otorhinolaryngol.* 2020;135:110136.

11. Yawn RJ, Dedmon MM, O'Connel BP, et al. Tympanic membrane perforation repair using porcine small intestinal submucosal grafting. *Otol Neurotol.* 2018;39(5):e332-e335.

12. Dontu P, Shaigany K, Eisenman DJ. Anatomic and audiometric outcomes of porcine intestinal submucosa compared to autologous fascia for tympanic membrane repair. Poster presented at: Combined Otolaryngology Spring Meetings, COSM 2022; April 27 – May 1, 2022; Dallas, TX.

® indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.



Biodesign®
OTOLOGIC BUTTERFLY GRAFT

For more detailed product information, including the Instructions for Use, which contains the indication statement, contraindications, precautions, and potential complications, visit www.everismedical.com.

