



Introducing the 6th generation Shaw Scalpel®

The SG6 is smarter, faster, and better than ever before!



- Faster calibration
- Improved blade responsiveness and ability to maintain set point temperature
- Ability to upgrade to the latest version of software via USB drive
- Video to support training and troubleshooting
- System memory for maintenance and diagnostics



Comparing the SG6 and the P8400



Characteristic	SG6	P8400
Handle and Blade Compatibility	Designed for 5800 series (single-use), made compatible with 5700 series	Designed for 5700 series, made compatible with 5800 series
Blade Temperature Control	Improved temperature control, reducing deviation range and improving performance	Up to 30°C variance from set point
Display	Visuals to convey blade temperature Secondary menu options for troubleshooting and in-service videos	
Software Upgrade	USB port to enable upgrades in the field	Not field upgradeable - requires return to C2Dx
Volume	Always available icon on the front of the controller	Adjusted by a radio-style knob on the rear of the controller
IV Pole Mounting	Both can be surface mounted by sitting on table or bench or secured to an IV pole via an integrated IV pole mount. Both systems utilize the same IV pole mount system.	
Temperature Range	Both have a temperature range from 70°C – 300°C	
Physical Size	Approximately 12 in. X 6 in. X 9 in. Weight: 12 lbs.	Approximately 8 in x 8 in x 12 in. Weight: 7.25 Lbs

A Better Experience for Everyone in the Room

- Quicker and smarter blade calibration and heating
- Less waiting
- See and experience less errors
- Greater product reliability
- Better visualization of blade status



Everis

555 East Eliza St, Ste. A
Schoolcraft, MI 49087 USA
t: 888 902 2239
www.everismedical.com

DOC-4323 Rev 1
Copyright © 2026 Everis

Everis™