


# Dilate salivary duct strictures.



## Advance<sup>®</sup>

### SALIVARY DUCT BALLOON CATHETER SET

Illustration by Lisa Clark

#### High-strength

The high-strength, minimally compliant material allows for the complete dilation of the stricture.

#### Minimally invasive

The flexible, low-profile design supports minimally invasive intervention through radial dilation.

#### Multi-modal visualization

The catheter is compatible with direct endoscopic visualization or ultrasound imaging based on physician preference.



# Advance®

## SALIVARY DUCT BALLOON CATHETER SET

Used for the dilation of strictures in the submandibular and parotid salivary ducts.

The balloon is constructed of an extra-thin wall and high-strength material that is intended to be used with sialendoscopic visualization in the submandibular salivary ducts.

### Components

- Advance Salivary Duct Balloon Catheter
- Wire guide

Order Number	Reference Part Number	Wire Guide Diameter inch	Inflated Balloon Diameter mm	Balloon Length cm	Preinflated Balloon Diameter mm
G34991	SDBS-14-30-1.5-2	.014	1.5	2	1
G34992	SDBS-14-30-2.0-2	.014	2.0	2	1
G34993	SDBS-14-30-2.5-2	.014	2.5	2	1
G34994	SDBS-14-30-3.0-2	.014	3.0	2	1

## Related products



### Sphere® Inflation Device

Order Number	Reference Part Number	Fluid Capacity mL	Maximum Pressure atm
G34903	EID-60-15	60	15

Some products or part numbers may not be available in all markets. Contact your local Everis sales representative or Customer Service for details.

### Everis

555 East Eliza St, Ste. A  
Schoolcraft, MI 49087 USA  
t: 888 902 2239  
www.everismedical.com

DOC-6573 Rev 2  
Copyright © 2026 Everis

# Everis™