


Everis™ Efer Bronchoscope

Size Chart




Bronchial Tubes for Adults



Part Number	Color	O.D. / I.D.	Length
BT-2000-3	Orange	13.2mm / 12.2mm	360mm / 36cm
BT-2101-3	Black	12.0mm / 11.0mm	360mm / 36cm
BT-2103-3	Red	10.0mm / 9.2mm	360mm / 36cm
BT-2105-3	Green	8.0mm / 7.0mm	360mm / 36cm
BT-2106-3	Blue	7.0mm / 6.5mm	360mm / 36cm

Tracheal Tubes for Adults



Part Number	Color	O.D. / I.D.	Length
BT-2210-3	Yellow	13.2mm / 12.2mm	260mm / 26cm
BT-2201-3	Black	12.0mm / 11.0mm	260mm / 26cm
BT-2203-3	Red	10.0mm / 9.2mm	260mm / 26cm
BT-2205-3	Green	8.0mm / 7.0mm	260mm / 26cm