



Everis™ Stent Placement Kit for Adults

The easy-to-use system places stents 8mm to 20mm in diameter and is equipped with:

- 4 stent placement sets color-coded identically with the bronchoscope tubes associated to them.
- 4 telescopic stent loaders for straight and Y stents.
- 1 double-action loader facilitating placement of both large straight and Y stents.
- A distal safety stop designed to prevent accidental damage to the trachea or bronchi during stent introducer tube insertions.
- A vision channel providing visual observation by an endoscope of the stent placement procedure.
- Ergonomic handles.



Stent Placement System



Ordering Information

Part Number	Description	Stent Diameter	Color
BR-2103-3	Stent Introducer for Adult	8mm-10mm	Red
BP-2103-3	Stent Pusher for Adult	8mm-10mm	Red
BI-2103-3	Stent Loader for Adult	8mm-10mm	Red
BR-2101-3	Stent Introducer for Adult	10mm-14mm	Black
BP-2101-3	Stent Pusher for Adult	10mm-14mm	Black
BI-2101-3	Stent Loader for Adult	10mm-14mm	Black
BR-2200-3	Stent Introducer for Adult	14mm-20mm	Yellow
BP-2200-3	Stent Pusher for Adult	14mm-20mm	Yellow
BI-2200-3	Stent Loader for Adult	14mm-20mm	Yellow
BR-2000-3	Stent Introducer for Adult	14mm-20mm and Y	Orange
BP-2000-3	Stent Pusher for Adult	14mm-20mm and Y	Orange
BI-2000-3	Stent Loader for Adult	14mm-20mm and Y	Orange
BI-2001-3	Double Action Stent Loader	Y and 14mm-20mm	Yellow/Orange

Learn more at www.everismedical.com

To place an order, email sales@everismedical.com

Please refer to the Intended Use and Instruction Manual for detailed information on insertion technique, sterilization, and postoperative care.

© indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.