



Everis™ Efer Bronchoscope for Adults

The bronchoscope is composed of:

- 5 color-coded ventilated bronchial tubes diameter 13.2mm / 12mm / 10mm / 8mm / 7mm.
- 4 color-coded tracheal tubes diameter 13.2mm / 12mm / 10mm / 8mm.
- 1 easily detachable universal base, secured locking by ring rotation.
- 1 compact base (allows for Tube-in-Tube Technique).
- 2 models of lateral obturators.
- Silicone caps.





Learn more at www.everismedical.com

To place an order, email sales@everismedical.com



The proximal end of all tubes fits directly into the ventilation connector when the base is disconnected.

Ordering Information

	Part Number	Description	Color
1	BB-2401-3	Lateral Port Adaptor with One Channel	
	BB-2402-3	Lateral Port Adaptor with Two Channels	
	BS-2101-3	Silicone Cap 12mm without Hole	Blue
	BS-2102-3	Silicone Cap 12mm for Instrument 2mm	Blue
2	BD-2410-3	Universal Base for Adults	
	BD-2501-3	Compact Base	
	BS-2300-3	Silicone Cap 25mm without Hole	Red
	BS-2320	Silicone Cap for Jet Ventilation Port	

	Part Number	Description		Color
3	BT-2000-3	Bronchial Tube for Adults	13.2mm x 12.2mm / 360mm	Orange
	BT-2101-3	Bronchial Tube for Adults	12mm x 11mm / 360mm	Black
	BT-2103-3	Bronchial Tube for Adults	10mm x 9.2mm / 360mm	Red
	BT-2105-3	Bronchial Tube for Adults	8mm x 7mm / 360mm	Green
	BT-2106-3	Bronchial Tube for Adults	7mm x 6.5mm / 360mm	Blue
	BT-2210-3	Tracheal Tube for Adults	13.2mm x 12.2mm / 260mm	Yellow
	BT-22010-3	Tracheal Tube for Adults	8mm x 7mm / 260mm	Green
	BT-2203-3	Tracheal Tube for Adults	10mm x 9.2mm / 260mm	Red
	BT-2205-3	Tracheal Tube for Adults	8mm x 7mm / 260mm	Green

Please refer to the Intended Use and Instruction Manual for detailed information on insertion technique, sterilization, and postoperative care.

® indicates U.S. trademark registration. All trademarks and/or images are the property of their respective owners or holders.