

Efer Bronchoscope Forceps

BIOPSY FORCEPS

2 mobile jaws & irrigation adapter

BPS-2001



BPS-2001

CROCODILE FORCEPS

2 mobile jaws & irrigation adapter

BPS-2002



BPS-2002

STENT REMOVAL FORCEPS

1 mobile jaw & irrigation adapter

BPS-2003



BPS-2003

HOOK SCISSORS FORCEPS

1 mobile jaw & irrigation adapter

BPS-2004



BPS-2004

Y STENT PLACEMENT FORCEPS

2 mobile jaws & irrigation adapter

BPS-2005



BPS-2005

STENT MEASUREMENT FORCEPS

2 mobile jaws & irrigation adapter

BPS-2006



BPS-2006

FENESTRATED GRASPING FORCEPS

1 fenestrated mobile jaw & irrigation adapter

BPS-2007



BPS-2007

REVERSE GRASPING FORCEPS

1 reverse mobile jaw & irrigation adapter

BPS-2008



BPS-2008

OPTICAL BIOPSY FORCEPS

2 mobile jaws

BPO-2000



OPTICAL CROCODILE FORCEPS

2 mobile jaws

BPO-2001



OPTICAL STENT REMOVAL FORCEPS

1 mobile jaw

BPO-2002



Ordering Information

Part Number	Description	Shaft Diameter	Head Diameter	Working Length
BPS-2001	Biopsy Forceps	1.5mm	3mm	60cm
BPS-2002	Crocodile Forceps	1.5mm	3mm	60cm
BPS-2003	Stent Removal Forceps	2.5mm	3.5mm	60cm
BPS-2004	Hook Scissors Forceps	1.5mm	3mm	60cm
BPS-2005	Y Stent Placement Forceps	5mm	5mm	60cm
BPS-2006	Stent Measurement Forceps		3mm	60cm
BPS-2007	Fenestrated Grasping Forceps	1.5mm	3mm	60cm
BPS-2008	Reverse Grasping Forceps	2.5mm	4.5mm	60cm
BPO-2000	Optical Biopsy Forceps		4mm	47cm
BPO-2001	Optical Crocodile Forceps		3mm	47cm
BPO-2002	Optical Stent Removal Forceps		4mm	47cm

Learn more at www.everismedical.com

To place an order, email sales@everismedical.com

Please refer to the Intended Use and Instruction Manual for detailed information on insertion technique, sterilization, and postoperative care.

Everis

555 East Eliza St, Ste. A
Schoolcraft, MI 49087 USA
t: 888 902 2239
www.everismedical.com

DOC-17164 Rev 1
Copyright © 2026 Everis

Everis[™]